

## N-Acetyl-5-hydroxytryptamine

Catalog No: #S1354


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	N-Acetyl-5-hydroxytryptamine
Brief Description	Inhibitors
Purification	99.00%
Target Name	5-HT Receptor inhibitor
Calculated MW	218.25
Formulation	C <sub>12</sub> H <sub>14</sub> N <sub>2</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: CC(=O)NCCC1=CNC2=C1C=C(O)C=C2 Pathways: GPCR/G Protein; Neuroscience Receptor: Sert (Sodium-dependent) Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Jang SW, et al. Proc Natl Acad Sci U S A. 2010 Feb 23; 107(8):3876-81.

Note: This product is for in vitro research use only