

Methyl 3,5-dihydroxybenzoate

Catalog No: #S1355


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Methyl 3,5-dihydroxybenzoate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	168.15
Formulation	C8H8O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:COC(=O)C1=CC(O)=CC(O)=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:Off-white Color PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lim Y, et al. Arch Pharm Res. 2011 Jul;34(7):1191-9.

Note: This product is for in vitro research use only