

TriMetazidine EP Impurity-D

Catalog No: #S1357


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TriMetazidine EP Impurity-D
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	198.22
Formulation	C10H14O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: COC1=C(OC)C(OC)=C(CO)C=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Belal TS, et al. J AOAC Int. 2014 Nov-Dec;97(6):1514-8.

Note: This product is for in vitro research use only