

2,3-Dichloroanisole


Catalog No: #S1358

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2,3-Dichloroanisole
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	177.03
Formulation	C7H6Cl2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=C(Cl)C(Cl)=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ramu L, et al. Magn Reson Chem. 2010 Aug;48(8):593-9.

Note: This product is for in vitro research use only