

## 2,3-Dihydrobenzofuranyl-5-acetic acid

Catalog No: #S1359


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	2,3-Dihydrobenzofuranyl-5-acetic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	178.18
Formulation	C10H10O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:OC(=O)CC1=CC2=C(OCC2)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Casalini C, et al. J Pharm Sci. 1980 Feb;69(2):164-7.

Note: This product is for in vitro research use only