

4-Nitro-3-trifluoromethyl aniline

Catalog No: #S1360


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4-Nitro-3-trifluoromethyl aniline
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	206.12
Formulation	C7H5F3N2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: NC1=CC(=C(C=C1)[N+](=O)C(F)(F)F
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Needle crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Saravanan S, et al. Spectrochim Acta A Mol Biomol Spectrosc. 2014; 121:685-97.

Note: This product is for in vitro research use only