

6-Hydroxypurine

Catalog No: #S1362


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	6-Hydroxypurine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	136.11
Formulation	C5H4N4O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:O=C1N=CNC2=C1NC=N2Pathways:OthersReceptor:Purine nucleoside phosphorylaseBoiling pt:Melting pt:Solubility:DMSO: 0.4 mg/mL (2.93 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Overington JP, et al. Nat Rev Drug Discov. 2006 Dec;5(12):993-6.

Note: This product is for in vitro research use only