

Cefalexine EP Impurity B

Catalog No: #S1364


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cefalexine EP Impurity B
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	214.24
Formulation	C ₈ H ₁₀ N ₂ O ₃ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CC1=C(N2[C@@H]([C@@H](C2=O)N)SC1)C(=O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Limited solubility Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Robin J, et al. Biotechnol Bioeng. 2003 Aug 5;83(3):353-60.

Note: This product is for in vitro research use only