

## Thioisonicotinamide


Catalog No: #S1365

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Thioisonicotinamide
Brief Description	Inhibitors
Purification	98%+
Target Name	Antibiotic
Calculated MW	138.19
Formulation	C6H6N2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: NC(=S)C1=CC=NC=C1 Pathways: Microbiology & virology Receptor: Tubercle bacillus Boiling pt: Melting pt: Solubility: Limited solubility Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Yamamoto K, et al. Saishin Igaku. 1967 Sep 10; 22(9): 2087-97.

Note: This product is for in vitro research use only