

3-Pyridylacetic acid hydrochloride

Catalog No: #S1366


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	3-Pyridylacetic acid hydrochloride
Brief Description	Inhibitors
Purification	98%+
Target Name	DPP-4 inhibitor
Calculated MW	173.6
Formulation	C7H8ClNO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: Cl.CC(=O)CC1=CN=CC=C1
Pathways: Proteases/Proteasome
Receptor: DPP-4
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Miyamoto Y, et al. J Med Chem. 2011 Feb 10;54(3):831-50.

Note: This product is for in vitro research use only