

Demeton-s-methyl sulfone

Catalog No: #S1367


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Demeton-s-methyl sulfone
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	262.28
Formulation	C6H15O5PS2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CCS(=O)(=O)CCSP(=O)(OC)OCPPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in WaterAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Burge WD. Appl Microbiol. 1969 Apr;17(4):545-50.

Note: This product is for in vitro research use only