

Methyl protocatechuate

Catalog No: #S1368


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Methyl protocatechuate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Adenosine Receptor; Apoptosis inducer
Calculated MW	168.15
Formulation	C ₈ H ₈ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC(=O)C1=CC(O)=C(O)C=C1
Pathways: Apoptosis; GPCR/G Protein Receptor: A2A; H₂O₂ induced apoptosis
Boiling pt: Melting pt: Solubility: Soluble in DMSO
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Zhang Z, et al. Neuroreport. 2015 Apr 15;26(6):367-73.

Note: This product is for in vitro research use only