

1H-1,2,4-Triazole-3-carboxylic acid


Catalog No: #S1371

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	1H-1,2,4-Triazole-3-carboxylic acid
Brief Description	Inhibitors
Purification	99.12%
Target Name	Others inhibitor
Calculated MW	113.07
Formulation	C3H3N3O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:OC(=O)C1=NC=NN1Pathways:OthersReceptor:isoprenyl transferasesBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. patent: US6221865

Note: This product is for in vitro research use only