

## L-5-Hydroxytryptophan

Catalog No: #S1379


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	L-5-Hydroxytryptophan
Brief Description	Inhibitors
Purification	99.00%
Target Name	Amino Acids and Derivatives inhibitor
Calculated MW	220.22
Formulation	C <sub>11</sub> H <sub>12</sub> N <sub>2</sub> O <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
SMILES: C1=CC2=C(C=C1O)C(=CN2)C[C@@H](C(=O)O)N  
Pathways: Others  
Receptor: Amino Acids and Derivatives  
Boiling pt: Melting pt: Solubility: Soluble in Water  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Overington JP, et al. Nat Rev Drug Discov. 2006 Dec;5(12):993-6.

Note: This product is for in vitro research use only