

Alodan

Catalog No: #S1380


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alodan
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	397.77
Formulation	C ₉ H ₆ Cl ₈
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: C1C(C(Cl)C(Cl)C(Cl)=C(Cl)C1(Cl)C2(Cl)C(Cl)C(Cl)C2)Cl
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. REHM WF, GARMS R. Z Tropenmed Parasitol. 1959 Aug;10:154-9.

Note: This product is for in vitro research use only