

## Biscomate


Catalog No: #S1383

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Biscomate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	192.26
Formulation	C <sub>8</sub> H <sub>4</sub> N <sub>2</sub> S <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: S=C=NC1=CC=C(C=C1)N=C=S  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: 129-131°C  
Solubility: Soluble in Ethanol  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Chodera L, et al. Pol Tyg Lek. 1973 Feb 26;28(9):338-40.

Note: This product is for in vitro research use only