

eIF4E Polyclonal Antibody

Catalog No: #40871

Package Size: #40871-1 50ul #40871-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	eIF4E Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	eIF4E Polyclonal Antibody detects endogenous levels of eIF4E protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human eIF4E around the non-phosphorylation site of S209.
Target Name	eIF4E
Other Names	EIF4E; EIF4EL1; EIF4F; Eukaryotic translation initiation factor 4E; eIF-4E; eIF4E; eIF-4F 25 kDa subunit; mRNA cap-binding protein
Accession No.	Swiss-Prot: P06730NCBI Gene ID: 1977
Uniprot	P06730
GeneID	1977;
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of COLO205 cells using eIF4E Polyclonal Antibody



Note: This product is for in vitro research use only