

Panthenol


Catalog No: #S1390

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Panthenol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Vitamin
Calculated MW	205.26
Formulation	C ₉ H ₁₉ NO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CC(C)(CO)C(O)C(=O)NCCCOP Pathways: Metabolism Receptor: VB5 Boling pt: 483.6°C Melting pt: 66-69°C Solubility: Free soluble in Water. Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Bissett DL. Clin Dermatol. 2009 Sep-Oct;27(5):435-45.

Note: This product is for in vitro research use only