

## Fluorescein

Catalog No: #S1392


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Fluorescein
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	332.32
Formulation	C20H12O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: OC1=CC2=C(C=C1)C1(OC(=O)C3=CC=CC=C13)C1=C(O2)C=C(O)C=C1 Pathways: Others Receptor: Ig kappa chain V-II region RPMI Boling pt: 620.8 Melting pt: 320C Solubility: Soluble in Water Appearance: Brownness Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Estevez-Torres A, et al. Anal Chem. 2007 Sep 25;.

Note: This product is for in vitro research use only