

Diatrizoic Acid

Catalog No: #S1395


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Diatrizoic Acid
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	613.92
Formulation	C11H9I3N2O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:CC(=O)NC1=C(I)C(C(O)=O)=C(I)C(NC(C)=O)=C1IPathways:OthersReceptor:OthersBoiling pt:614.1CMelting pt:>300CSolubility:Soluble in WaterAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Rastogi T, et al. Sci Total Environ. 2014 Jun 1;482-483:378-88.

Note: This product is for in vitro research use only