

Agipnon

Catalog No: #S1396


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Agipnon
Brief Description	Inhibitors
Purification	99.00%
Target Name	GABA Receptor antagonist
Calculated MW	155.19
Formulation	C8H13NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CCC1(C)CC(=O)NC(=O)C1Pathways:Membrane transporter/Ion channel;

NeuroscienceReceptor:GABARBoiling pt:282C?Melting pt:126-127CSolubility:Ethanol: 31 mg/mL (199.75 mM); DMSO: 31 mg/mL (199.75

mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Mistry DK, et al. Exp Physiol, 1990, 75(2),199-209.

Note: This product is for in vitro research use only