

Sodium montmorillonite

Catalog No: #S1397

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sodium montmorillonite
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	282.21
Formulation	Al ₂ O ₉ Si ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Metabolism

system SMILES: O.[O-].[O-].[O-].[Al+3].[Al+3].O=[Si]=O.O=[Si]=O.O=[Si]=O.O=[Si]=O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Khaki Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Joshi PC, Aldersley MF. J Mol Evol. 2013 Jun;76(6):371-9

Note: This product is for in vitro research use only