

Climbazole

Catalog No: #S1403


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Climbazole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	292.77
Formulation	C15H17ClN2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CC(C)(C)C(=O)C(OC1=CC=C(Cl)C=C1)N1C=CN=C1
 Pathways: Others
 Receptor: Others
 Boiling pt: 447.5°C
 Melting pt: 96-100°C
 Solubility: Ethanol: 59 mg/mL (201.53 mM); DMSO: 59 mg/mL (201.53 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Prez-Rivera AA, et al. Mutat Res. 2009 Jan;672(1):27-39.

Note: This product is for in vitro research use only