

Cefditoren pivoxil

Catalog No: #S1405


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cefditoren pivoxil
Brief Description	Inhibitors
Purification	93.80%
Target Name	Antibacterial
Calculated MW	620.72
Formulation	C25H28N6O7S3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1=C(SC=N1)/C=C\C2=C(N3[C@@H]([C@@H](C3=O)NC(=O)/C(=NOC)/C4=CSC(=N4)N)SC2)C(=O)OCOC(=O)C(C)C)C Pathways: Microbiology & virology Receptor: PBPs Boiling pt: Melting pt: 207-209°C Solubility: DMSO: 100 mg/mL (161.1 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Martin SI, Kaye KM. Infect Dis Clin North Am. 2004 Sep;18(3):603-19, ix.

Note: This product is for in vitro research use only