

Pentetrazol


Catalog No: #S1413

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pentetrazol
Brief Description	Inhibitors
Purification	98.00%
Target Name	GABA Receptor antagonist
Calculated MW	138.17
Formulation	C6H10N4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: C1CCN2N=NN=C2CC1
Pathways: Membrane transporter/Ion channel; Neuroscience
Receptor: GABA
Boiling pt: 194°C
Melting pt: 59-61°C
Solubility: DMSO: > 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Squires RF, et al. Life Sci. 1984, 35 (14): 1439C44.

Note: This product is for in vitro research use only