

Mesna

Catalog No: #S1414


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mesna
Brief Description	Inhibitors
Purification	99.93%
Target Name	Others
Calculated MW	164.18
Formulation	C ₂ H ₅ NaO ₃ S ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Other SMILES:[Na+].[O-]S(=O)(=O)CC Pathways:Others Receptor:Others Boiling pt:Melting pt:>240C? Solubility:Water: 7 mg/mL (42.63 mM); DMSO: 33 mg/mL (200.99 mM) Appearance:White Powder Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kalcioğlu MT, et al. Am J Otolaryngol. 2014 May-Jun;35(3):357-61.

Note: This product is for in vitro research use only