

Fosfosal


Catalog No: #S1427

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fosfosal
Brief Description	Inhibitors
Purification	97.00%
Target Name	COX inhibitor
Calculated MW	218.1
Formulation	C7H7O6P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: OC(=O)C1=CC=CC=C1OP(=O)(O)O
Pathways: Immunology/Inflammation; Neuroscience
Receptor: COX
Boiling pt: 482.3°C
Melting pt: 162-166°C
Solubility: Water: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Gmez LA, et al. Eur J Pharmacol. 1988 Mar 15;147(3):459-63.

Note: This product is for in vitro research use only