

## Betamipron

Catalog No: #S1430


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                   |   |
|-------------------|---|
| Product Name      | Betamipron                                      |
| Brief Description | Inhibitors                                      |
| Purification      | 99.00%  |
| Target Name       | Others inhibitor                                |
| Calculated MW     | 193.2   |
| Formulation       | C10H11NO3                                       |
| Storage           | 3 years -20°C powder; 2 years -80°C in solvent; |

## Images



## Product Description

Research Area: Immunology  
 SMILES: OC(=O)CCNC(=O)C1=CC=CC=C1  
 Pathways: Others Receptor: Others  
 Boiling pt: 469.2°C  
 Melting pt: 132°C  
 Solubility: Ethanol: 39 mg/mL (201.86 mM); DMSO: 39 mg/mL (201.86 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Ohashi N, et al. J Infect ChemOthers, 2005, 11(1), 24-31.

Note: This product is for in vitro research use only