

Khellin

Catalog No: #S1431


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Khellin
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	260.25
Formulation	C ₁₄ H ₁₂ O ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
SMILES: COC1=C2OC(C)=CC(=O)C2=C(OC)C2=C1OC=C2
Pathways: Others
Receptor: Others
Boiling pt: 482.1°C
Melting pt: 150-154°C
Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Vrzal R, et al. PLoS One. 2013 Sep 19;8(9):e74917.

Note: This product is for in vitro research use only