

Dexpanthenol

Catalog No: #S1437


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dexpanthenol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Vitamin inhibitor
Calculated MW	205.25
Formulation	C9H19NO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:CC(C)(CO)[C@H](C(=O)NCCCO)OPathways:MetabolismReceptor:VB5Boiling pt:118-120CMelting pt:Solubility:Soluble in WaterAppearance:Transparent Paste LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ebner F, et al. Am J Clin Dermatol. 2002; 3(6):427-33.

Note: This product is for in vitro research use only