

## Butoconazole nitrate


Catalog No: #S1438

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Butoconazole nitrate
Brief Description	Inhibitors
Purification	99.50%
Target Name	Estrogen/progestogen Receptor inhibitor
Calculated MW	474.79
Formulation	C19H18Cl3N3O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: C1C=CC=C(CCC(CN2C=CN=C2)SC2=C(Cl)C=CC=C2Cl)C=C1 Pathways: Endocrinology/Hormones Receptor: ER Boling pt: 566.9C Melting pt: 159C Solubility: DMSO: 95 mg/mL (200.08 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Pfaller MA, et al. *Diagn Microbiol Infect Dis*. 1990 Jan-Feb;13(1):31-5.

Note: This product is for in vitro research use only