

Ethacridine lactate monohydrate

Catalog No: #S1440


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ethacridine lactate monohydrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	361.4
Formulation	C ₁₈ H ₂₃ N ₃ O ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: O.CC(O)C(O)=O.CCOC1=CC2=C(N)C3=C(C=C(N)C=C3)N=C2C=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Water: 5 mg/mL (13.83 mM); DMSO: 72 mg/mL (199.23 mM) Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Junka A, et al. Int Wound J, 2013.

Note: This product is for in vitro research use only