

Pizotifen

Catalog No: #S1442


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pizotifen
Brief Description	Inhibitors
Purification	99.00%
Target Name	5-HT Receptor antagonist
Calculated MW	295.44
Formulation	C19H21NS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CN1CCC(CC1)=C1C2=C(CCC3=CC=CC=C13)SC=C2
Pathways: GPCR/G Protein;
Neuroscience
Receptor: 5-HT
Boiling pt: Melting pt: 140-142°C
Solubility: Soluble in DMSO
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Rahman SM, et al. Pak J Pharm Sci. 2010 Oct;23(4):435-41.

Note: This product is for in vitro research use only