

## Dasatinib

Catalog No: #S1448


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Dasatinib
Brief Description	Inhibitors
Purification	99.62%
Target Name	Bcr-Abl inhibitor; c-Kit inhibitor; Ephrin Receptor inhibitor; Src inhibitor
Calculated MW	488.01
Formulation	C22H26ClN7O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer  
 SMILES: CC1=NC(=CC(NC2=NC=C(S2)C(=O)NC2=C(C)C=CC=C2Cl)=N1)N1CCN(CCO)CC1  
 Pathways: Angiogenesis; Tyrosine Kinase/Adaptors; Cytoskeletal Signaling  
 Receptor: Abl; c-Kit (D816V); c-Kit (wt); EphA2; Src  
 Boling pt: Melting pt: 275-286°C  
 Solubility: DMSO: 98 mg/mL (200.81 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

- O'Hare T, et al. Y Res. 2005, 65(11), 4500-4505.

Note: This product is for in vitro research use only