

Toremifene citrate

Catalog No: #S1464


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Toremifene citrate
Brief Description	Inhibitors
Purification	99.50%
Target Name	Estrogen/progestogen Receptor modulator
Calculated MW	598.08
Formulation	C32H36ClNO8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer
 SMILES: CN(C)CCOC1=CC=C(C=C1)/C(=C(/CCCl)\C2=CC=CC=C2)/C3=CC=CC=C3.C(C(=O)O)C(CC(=O)O)(C(=O)O)O
 Pathways: Endocrinology/Hormones
 Receptor: ER
 Boiling pt: 146-148°C
 Melting pt: 160-162°C
 Solubility: DMSO: 100 mg/mL (167.2 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Locci P, et al. J Interferon Cytokine Res. 2001 Nov;21(11):961-70.

Note: This product is for in vitro research use only