

Nadifloxacin

Catalog No: #S1466


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nadifloxacin
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	360.39
Formulation	C ₁₉ H ₂₁ FN ₂ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CC1CCC2=C3N1C=C(C(O)=O)C(=O)C3=CC(F)=C2N1CCC(O)CC1
 Pathways: Others
 Receptor: Others
 Boiling pt: 624.9°C
 Melting pt: 245-247°C
 Solubility: DMSO: 20 mg/mL (55.49 mM)
 Appearance: Off-white Color Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Jacobs MR, et al. Expert Opin PharmacOthers, 2006, 7(14), 1957-1966.

Note: This product is for in vitro research use only