

Levocetirizine

Catalog No: #S1469


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Levocetirizine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Histamine Receptor antagonist
Calculated MW	388.89
Formulation	C21H25ClN2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:C1CN(CCN1CCOCC(=O)O)[C@H](C2=CC=CC=C2)C3=CC=C(C=C3)ClPathways:GPCR/G Protein; Immunology/Inflammation; NeuroscienceReceptor:H1 receptorBoiling pt:542.1CMelting pt:205-208CSolubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Pasquali M, et al. Clin Exp Allergy. 2006 Sep;36(9):1161-7.

Note: This product is for in vitro research use only