

Fexofenadine hydrochloride

Catalog No: #S1470


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fexofenadine hydrochloride
Brief Description	Inhibitors
Purification	99.50%
Target Name	Histamine Receptor antagonist
Calculated MW	538.12
Formulation	C32H40ClNO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:[H+].[Cl-].CC(C)(C(O)=O)C1=CC=C(C=C1)C(O)CCCN1CCC(CC1)C(O)(C1=CC=CC=C1)C1=CC=CC=C1Pathways:GPCR/G Protein;

Immunology/Inflammation; NeuroscienceReceptor:H1 receptorBoiling pt:697.3CMelting pt:Solubility:DMSO: 10 mMAppearance:White

PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Abbas MN, et al. Anal Sci. 2004 Aug;20(8):1137-42.

Note: This product is for in vitro research use only