

Trilostane

Catalog No: #S1473


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trilostane
Brief Description	Inhibitors
Purification	99.50%
Target Name	Dehydrogenase inhibitor
Calculated MW	329.43
Formulation	C ₂₀ H ₂₇ NO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@@H]2O)CC[C@]45[C@@]3(CC(=C([C@H]4O5)O)C#N)CPathways:MetabolismReceptor:3 β-hydroxysteroid dehydrogenaseBoiling pt:516.5CMelting pt:264CSolubility:DMSO: 65 mg/mL (197.31 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Puddefoot JR, et al. Expert Opin PharmacOthers. 2006 Dec;7(17):2413-9.

Note: This product is for in vitro research use only