

Emp Polyclonal Antibody

Catalog No: #40880

Package Size: #40880-1 50ul #40880-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Emp Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Emp Polyclonal Antibody detects endogenous levels of Emp protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Emp.
Target Name	Emp
Other Names	MAEA; EMP; HLC10; PIG5; Macrophage erythroblast attacher; Cell proliferation-inducing gene 5 protein; Erythroblast macrophage protein; Human lung cancer oncogene 10 protein; HLC-10
Accession No.	Swiss-Prot: Q7L5Y9NCBI Gene ID: 10296
Uniprot	Q7L5Y9
GeneID	10296;
SDS-PAGE MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using Emp Polyclonal Antibody

Note: This product is for in vitro research use only