

Fudosteine

Catalog No: #S1481


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fudosteine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	179.23
Formulation	C6H13NO3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:C(CO)CSC[C@@H](C(=O)O)NPathways:OthersReceptor:MUC5AC mucinBoiling pt:Melting pt:200-202CSolubility:Water: 36 mg/mL (200.84 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Rhee CK, et al. Eur Respir J. 2008 Nov;32(5):1195-202.

Note: This product is for in vitro research use only