

Perindoprilat

Catalog No: #S1484


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Perindoprilat
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR inhibitor; MRP inhibitor; RAAS inhibitor
Calculated MW	340.41
Formulation	C17H28N2O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: CCCC[C@@H](C(=O)O)N[C@@H](C(=O)N1[C@H]2CCCC[C@H]2C[C@H]1C(=O)O Pathways: Endocrinology/Hormones;

Immunology/Inflammation; Neuroscience Receptor: AChR; MRP; ACEBoling pt: Melting pt: 153-155°C Solubility: Soluble in DMSO Appearance: White

Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Alfakih K, Hall AS. Expert Opin PharmacOther. 2006 Jan;7(1):63-71.

Note: This product is for in vitro research use only