

## Perindopril erbumine

Catalog No: #S1484L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Perindopril erbumine
Brief Description	Inhibitors
Purification	98.00%
Target Name	MRP inhibitor; RAAS inhibitor
Calculated MW	441.61
Formulation	C23H43N3O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

SystemSMILES: CCC[C@@H](C(=O)OCC)N[C@@H](C)C(=O)N1[C@H]2CCCC[C@H]2[C@H]1C(=O)O.CC(C)(C)N  
 Pathways: Endocrinology/Hormones; Immunology/Inflammation  
 Receptor: MRP1; ACE  
 Boiling pt: 537.4°C  
 Melting pt: 126-128°C  
 Solubility: Ethanol: 88 mg/mL (199.27 mM); Water: 88 mg/mL (199.27 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
 Stock solution can be stored below -20 °C for several month.

## References

1. Ceconi C, et al. Eur J Pharmacol, 2007, 577(1-3), 1-6.

Note: This product is for in vitro research use only