

Gestodene

Catalog No: #S1490


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gestodene
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	310.43
Formulation	C ₂₁ H ₂₆ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CC[C@]12CC[C@H]3[C@H]([C@@H]1C=C[C@]2(C#C)O)CCC4=CC(=O)CC[C@H]34Pathways:Endocrinology/HormonesReceptor:

ERBoiling pt:462.7C?Melting pt:190-192CSolubility:Ethanol: 8 mg/mL (25.77 mM); DMSO: 62 mg/mL (199.72 mM)Appearance:White

PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lidegaard ?, et al. BMJ. 2011 Oct 25;343:d6423.

Note: This product is for in vitro research use only