

## Alarelin

Catalog No: #S1494


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Alarelin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others agonist
Calculated MW	1167.32
Formulation	C56H78N16O12
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:CCNC(=O)[C@@H]1CCCN1C(=O)[C@H](CCC=C(N)N)NC(=O)[C@H](CC(C)C)NC(=O)[C@@H](C)NC(=O)[C@H](CC2=CC=C(C=C2)O)NC(=O)[C@H](CO)NC(=O)[C@H](CC3=CNC4=CC=CC=C43)NC(=O)[C@H](CC5=CN=CN5)NC(=O)[C@@H]6CCC(=O)N6.CC(=O)O

Pathways :OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Suocheng Weia, et al. Small Ruminant Research.2012,105 (1-3): 237-243.

Note: This product is for in vitro research use only