

Liranaftate

Catalog No: #S1497


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Liranaftate
Brief Description	Inhibitors
Purification	99.82%
Target Name	Antifungal inhibitor
Calculated MW	328.43
Formulation	C18H20N2O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: COC1=CC=CC(=N1)N(C)C(=S)OC1=CC2=C(CCCC2)C=C1
 Pathways: Microbiology & virology
 Receptor: squalene epoxidase
 Boiling pt: 462.5°C
 Melting pt: 98.5-99.5°C
 Solubility: DMSO: >5mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Maruyama N, et al. Med Mycol J. 2012;53(2):129-33.

Note: This product is for in vitro research use only