

Norethindrone

Catalog No: #S1514


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Norethindrone
Brief Description	Inhibitors
Purification	99.33%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	298.42
Formulation	C ₂₀ H ₂₆ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@]2(C#C)O)CCC4=CC(=O)CC[C@H]34Pathways:Endocrinology/HormonesReceptor:ER
 Boling pt:447C?Melting pt:205-206CSolubility:Ethanol: 5 mg/mL (16.75 mM); DMSO: 60 mg/mL (201.05 mM)Appearance:White PowderRemark:For
 obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C
 for several month.

References

1. Mendoza-Rodriguez CA, et al. Contraception. 1999 May;59(5):339-43.

Note: This product is for in vitro research use only