

Bromindione


Catalog No: #S1515

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bromindione
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	301.14
Formulation	C15H9BrO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: BrC1=CC=C(C=C1)C1C(=O)C2=CC=CC=C2C1=O
Pathways: Others
Receptor: Others
Boiling pt: 450.1°C
Melting pt: 139-142°C
Solubility: Soluble in DMSO
Appearance: Purple Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

1. SINGER MM, et al. JAMA. 1962 Jan 13;179:150-2.

Note: This product is for in vitro research use only