

## Benicar

Catalog No: #S1518


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Benicar
Brief Description	Inhibitors
Purification	99.50%
Target Name	RAAS inhibitor
Calculated MW	558.59
Formulation	C29H30N6O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

System SMILES: CCCC1=NC(=C(N1CC1=CC=C(C=C1)C1=CC=CC=C1C1=NNN=N1)C(=O)OCC1=C(C)OC(=O)O1)C(C)(C)O Pathways: Endocrinology/Hormones Receptor: angiotensin II receptor Boiling pt: 804.2°C Melting pt: 180°C Solubility: DMSO: 89 mg/mL (159.32 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Koike H, et al. J Hypertens Suppl. 2001 Jun;19(1):S3-14.

Note: This product is for in vitro research use only