

## Iopromide

Catalog No: #S1522


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Iopromide
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	791.12
Formulation	C18H24I3N3O8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Other SMILES: COCC(=O)NC1=C(I)C(C(=O)NCC(O)CO)=C(I)C(C(=O)N(C)CC(O)CO)=C1 I Pathways: Others Receptor: Others Boiling pt: 840.9°C Melting pt: Solubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Palkowitsch PK, et al. Acta Radiol. 2014 Jul;55(6):707-14.

Note: This product is for in vitro research use only